

# Profiling and Relative Quantitation of Phosphoinositides by Multiple Precursor Ion Scanning Based on Phosphate Methylation and Isotopic Labeling

Tanxi Cai,<sup>†,‡,⊥</sup> Qingbo Shu,<sup>†,‡,⊥</sup> JunJie Hou,<sup>§</sup> Peibin Liu,<sup>†</sup> Lili Niu,<sup>†</sup> Xiaojing Guo,<sup>†</sup> Charles C. Liu,<sup>||</sup> and Fuquan Yang<sup>\*,†</sup>

<sup>†</sup>Laboratory of Protein and Peptide Pharmaceuticals and Laboratory of Proteomics, Institute of Biophysics, Chinese Academy of Sciences, Beijing 100101, China


<sup>‡</sup>University of Chinese Academy of Sciences, Beijing 100049, China

<sup>§</sup>National Laboratory of Biomacromolecules, Institute of Biophysics, Chinese Academy of Sciences, Beijing 100101, China

<sup>||</sup>ASPEC Technologies Limited, Beijing, China

## Supporting Information

**ABSTRACT:** Phosphoinositides, the phosphorylated derivatives of phosphatidylinositol (PtdIns), are key regulators of many fundamental biological processes, including cell growth, proliferation, and motility. Here, we present a novel method for rapid, sensitive, and simultaneous profiling of phosphatidylinositol trisphosphate (PtdInsP<sub>3</sub>), phosphatidylinositol bisphosphate (PtdInsP<sub>2</sub>), and phosphatidylinositol phosphate (PtdInsP) of different fatty acid compositions. This method is based on a technique called “charged diacylglycerol fragment ion-specific multiple precursor ion scanning” (DAG<sup>+</sup>-specific MPIS), coupled with prior phosphate methylation. Using DAG<sup>+</sup>-specific MPIS, we were able to identify 32 PtdIns, 28 PtdInsP, 30 PtdInsP<sub>2</sub>, and 3 PtdInsP<sub>3</sub> molecular species from bovine brain extracts or prostatic cancer cell lines in an efficient and time-saving manner. Our analysis revealed a large range of fatty acyl compositions in phosphoinositides not obtained previously from mammalian samples. We also developed a method that involves isotopic labeling of endogenous phosphoinositides with deuterated diazomethane (CD<sub>2</sub>N<sub>2</sub>) for quantitation of phosphoinositides. CD<sub>2</sub>N<sub>2</sub> was generated in situ through acid-catalyzed H/D exchange and methanolysis of trimethylsilyl diazomethane (TMS-diazomethane). Phosphoinositides, extracted from a PC3 prostatic cancer cell line, were labeled either with CH<sub>2</sub>N<sub>2</sub> or CD<sub>2</sub>N<sub>2</sub> and mixed in known proportions for DAG<sup>+</sup>-specific MPIS-based mass spectrometry (MS) analysis. The results indicate that isotopic labeling is capable of providing accurate quantitation of PtdInsP<sub>3</sub>, PtdInsP<sub>2</sub>, and PtdInsP with adequate linearity as well as high reproducibility with an average coefficient variation of 18.9%. More importantly, this new methods excluded the need for multiple phosphoinositide internal standards. DAG<sup>+</sup>-specific MPIS and isotopic labeling based MS analysis of phosphoinositides offers unique advantages over existing approaches and presents a powerful tool for research of phosphoinositide metabolism.



Phosphoinositides, the polyphosphorylated derivatives of phosphatidylinositol (PtdIns), are found in virtually all cellular membranes, and they predominantly reside on the cytoplasmic side of the cell membrane. Reversible phosphorylation at the hydroxyl groups 3, 4, and 5 of the inositol ring in PtdIns by a number of different kinases generates seven distinct phosphoinositides, defined here as PtdIns(3)P, PtdIns(4)P, PtdIns(5)P, PtdIns(3,4)P<sub>2</sub>, PtdIns(3,5)P<sub>2</sub>, PtdIns(4,5)P<sub>2</sub>, and PtdIns(3,4,5)P<sub>3</sub>, contributing to their signaling diversity.<sup>1</sup> Though present in low abundance in mammalian membranes, phosphoinositides play pivotal role in numerous cellular processes, including cell growth, survival, and motility.<sup>2</sup> Aberrant phosphoinositide signaling has been implicated in numerous human diseases, such as cancer, neurological disorders, diabetes, and others.<sup>3,4</sup>

Despite the clear importance of phosphoinositides in biological signaling, there is still little known about the precise molecular species of phosphoinositides and the exact cellular amounts of these lipids across different cells, tissues, or organisms. However, measuring the phosphoinositides present a major challenge, as they are present at only very low levels inside the cell. For example, of the total cellular phospholipids content, PtdIns(4,5)P<sub>2</sub> and PtdIns(3,4,5)P<sub>3</sub>, are present at less than 1% and 0.05%, respectively.<sup>5</sup> In addition, receptor displacement assays and radiometric measurements of deacy-

**Received:** August 3, 2014

**Accepted:** November 26, 2014

**Published:** November 26, 2014



lated headgroup products of phosphoinositides are laborious and fail to elucidate the exact fatty acyl composition, in spite of the high sensitivity of such assays.<sup>6,7</sup>

Recent advances in separation technology and mass spectrometry have revolutionized the characterization and quantification of the phosphoinositides, allowing for measurements with high sensitivity as well as high-throughput setup, and enable analysis across a large dynamic range.<sup>8–10</sup> However, despite considerable progress in recent years, these powerful approaches have been hampered by major issues in terms of chemical stability and sensitivity associated with the measurement and quantification of phosphoinositides. Recently, Clark et al. used trimethylsilyl-diazomethane as a reagent that methylates phosphate groups of phosphoinositides from total cell or tissue extracts.<sup>11</sup> Such approach offers a novel route to solving the major problems mentioned above. This method was also shown to be suitable for determining phosphoinositides levels and fatty acid composition in small cell and tissue samples. The general outline of this method, however, is time-consuming, as it first has to run neutral loss scans (NLSs) of the methylated inositol phosphate headgroup, in order to identify the particular molecular species of phosphoinositide, so that phosphoinositides can be accurately quantified by multiple reaction monitoring (MRM)-based mass spectrometry analysis. In addition, such strategy is prone to biased losses of information on phosphoinositide molecular species of low abundance during NLSs analysis. This is really true in cases where peaks of phosphoinositide species in the NLS spectrum may be barely above the noise threshold. Furthermore, the comparative quantitation of phosphoinositides molecular species is usually achieved by the inclusion of internal standards (ISDs) and analyzed in separate runs. However, the stability of the electrospray might differ at times, and the potential difference in the ionization of phosphoinositides of different fatty acid composition and degree of unsaturation may prone to errors in quantification. Furthermore, for most phosphoinositide species, the relevant ISDs are not commercially available and therefore need to be synthesized chemically. The storage of ISDs in solution for extended periods of time is also problematic, especially for ISDs with double bonds due to susceptibility to rapid oxidation. To overcome limitations of such nature, a simple and straightforward approach would be to incorporate a given heavy nuclei (e.g., <sup>2</sup>H, <sup>13</sup>C, or <sup>15</sup>N) into endogenous phosphoinositides by metabolic or chemical labeling. The light and heavy populations of phosphoinositides are sufficiently equivalent chemically, but also remain distinguishable by mass spectrometry (MS). If two populations are mixed, the phosphoinositide abundances can be directly determined from the relative MS signal intensities. In support of our choice of labeling, numerous isotopic labeling procedures, such as SILAC (stable isotope labeling by amino acids in cell culture) or iTRAQ (isobaric tag for relative and absolute quantitation), have been established for the study of protein mixtures and are known for their robust and accurate quantification of protein expression levels.<sup>12,13</sup> In this study, we present an integrated method that combines charged diacylglycerol fragment ion-specific multiple precursor ion scanning (DAG<sup>+</sup>-specific MPIS), direct infusion mass spectrometry (DIMS), and isotopic labeling through phosphate methylation for rapid profiling and relative quantitation of phosphoinositides from various biological samples.


## ■ EXPERIMENTAL SECTION

**Chemicals and Reagents.** Trimethylsilyl diazomethane (IUPAC name, (diazomethyl)-trimethylsilane; abbreviated here as TMS-diazomethane) (2 M in hexanes), deuterium oxide (D<sub>2</sub>O), deuterium-labeled methanol (methanol-*d*<sub>4</sub>, MeOD), and deuterium chloride (DCl, 20 wt % in D<sub>2</sub>O) were purchased from Acros Organics. Acetic acid and hydrochloric acid (HCl) of analytical grade were purchased from Fisher Scientific. Chloroform (SpS grade, CHCl<sub>3</sub>), methanol (SpS grade, MeOH), water (ultragradient grade, H<sub>2</sub>O), and methyl *tert*-butyl ether (ultragradient grade, MTBE) were from J.T. Baker. The phosphoinositides sodium salt from bovine brain was from Sigma-Aldrich. Synthetic 17:0/C20:4-PtdIns(3,4,5)P<sub>3</sub>, 17:0/C20:4-PtdIns(4,5)P<sub>2</sub>, and 17:0/C20:4-PtdIns(4)P were from Avanti Polar Lipids. All the reagents used for cell culture were from Invitrogen unless stated otherwise. All other reagents were from Sigma-Aldrich unless specified.

**Cell Lines and Cell Cultures.** The prostate cancer cell lines PC3, Du145, and C4-2B were obtained from ATCC. The PC3, Du145 and C4-2B cells were cultured in DMEM medium with 10% (v/v) FBS and 1% (v/v) pen/strep. All the cells were maintained at 37 °C in a 5% CO<sub>2</sub> incubator.

**Lipid Extraction.** The lipid extraction was performed using 2 mL polypropylene safe-lock tubes. Two-step lipid extraction was performed as described previously with some modifications.<sup>14</sup> Briefly, cells were collected after digestion with trypsin and subjected to cell counting assays using the Neubauer plate. Prior to lipid extraction, ~5 × 10<sup>5</sup> cells were centrifuged (13 000g, 5 min, 4 °C) and the pellets were washed for two times with 1 mL of cold 5% TCA/1 mM EDTA. The supernatants were then discarded and neutral lipids were extracted by incubating with 1 mL of chloroform/methanol (1:2, v/v) and vortexing for 3–4 times over 10 min at room temperature (RT). After further centrifugation (13 000g, 2 min at 4 °C) and carefully removing the supernatant, 10 pmol of PtdIns(4)P 37:4, 10 pmol of PtdIns(4,5)P<sub>2</sub> 37:4, and 10 pmol of PtdIns(3,4,5)P<sub>2</sub> 37:4 were spiked into the remaining pellet. The pellet was resuspended in 1365 μL of MTBE/MeOH (MeOD, in case of isotopic labeling)/2 N HCl (DCl, in case of isotopic labeling) (200:60:13, v/v/v), incubated at RT for 15 min, and vortexed every 5 min for 30 s. An amount of 250 μL of 0.1 N HCl (or DCl) was added followed by 5 min of vortexing. Phase separation was induced by centrifugation (6500g, 2 min at 4 °C). The upper organic phase was immediately transferred to a new tube and added with 500 μL of prederivatization wash solution (MTBE/MeOH or MeOD/0.01 N HCl or DCl, v/v/v = 100:30:25, low phase). After a brief vortex and centrifugation (6500g, 2 min at 4 °C), the upper phase was transferred into another new tube for derivatization.

**Derivatization of Lipids.** TMS-diazomethane at 2 M in hexane (50 μL) was added into the lipid extracts prepared as described above (~1 mL of “upper phase”) and incubated for 20 min at RT. The reactions were terminated by adding 6 μL of glacial acetic acid, causing an easily visible color change. An amount of 500 μL of postderivatization wash solution (MTBE/MeOH/H<sub>2</sub>O, v/v/v = 100:30:25, low phase) was then added into the organic solution. The samples were then mixed and centrifuged (1500g, 3 min), and the upper phase was collected. The wash step was repeated once, and the final upper phase was then dried in SpeedVac and resuspend in 80 μL of



**Figure 1.** MS/MS analysis of completed, overmethylated 18:0/20:4-PtdInsP<sub>2</sub> and 17:0/20:4-PtdInsP<sub>3</sub>. (A and B) Collision-induced dissociation of the completed ( $m/z$  of 1117) or overmethylated ( $m/z$  of 1131) 18:0/20:4-PtdInsP<sub>2</sub> produce the same fragment ion at  $m/z$  of 627, which was corresponding to the charged 18:0/20:4-DAG. (C and D) Collision-induced dissociation of the completed ( $m/z$  of 1211) or overmethylated ( $m/z$  of 1225) 17:0/20:4-PtdInsP<sub>3</sub> produce the same fragment ion at  $m/z$  of 613, which was corresponding to the charged 17:0/20:4-DAG. Characteristic fragment ions and corresponding neutral losses are indicated.

chloroform/methanol (v/v) with 5 mM ammonium acetate before mass spectrometer analysis.

**Stable Isotopes Labeling on Lipids.** The labeling of lipids with stable isotopes was performed by replacing the MeOH, H<sub>2</sub>O, and HCl with MeOD, D<sub>2</sub>O, and DCl, respectively, during the lipid extraction and derivatization procedures.

**Mass Spectrometric Analysis.** Mass spectrometric analyses were performed on the TSQ Vantages instrument (ThermoFisher Scientific, Bremen), equipped with a robotic nanoflow ion source Triversa Nanomate (Advion Biosciences) using chips with spraying nozzles with a diameter of 5.5  $\mu$ m. The ion source was controlled using the ChipSoft Software (Advion Biosciences), with the following settings: ionization voltage was set to 1.25 kV, gas backpressure to 0.5 psi in the positive ion mode, 1  $\mu$ L air gap before chip and prewetting one time was enabled, and vent headspace was disabled. The temperature of the ion transfer capillary and the s-lens rf amplitude were 190  $^{\circ}$ C and 217. DAG<sup>+</sup>-specific MPIS was carried out at a CE of 40 eV, using a dwelling time of 500 ms at a step size of 0.2 Da at unit resolution of the Q1 quadrupole. For every precursor ion scanning, 24 spectra were acquired. The Q1 mass range was set to  $m/z$  700–1400. MS<sup>3</sup> analyses were performed on LTQ-Orbitrap XL instrument (ThermoFisher Scientific, Bremen). All MS<sup>3</sup> spectra were collected using normalized collision energy (setting, 35%) and an isolation window of 3  $m/z$ .

For DAG<sup>+</sup>-specific MPIS, we constructed a list of daughter ions that putatively correspond to charged diacylglycerols (DAG<sup>+</sup>) with common fatty acyl groups (see Supporting Information Table S1 for a list of putative DAG<sup>+</sup> and their associated masses). A scheduled DAG<sup>+</sup>-specific MPIS was then performed for a rapid screening of phosphoinositide molecular species from different biological sample in 10 min (Supporting Information Figure S1).

#### Quantitative Monitoring of Changes in Lipid Profile.

Relative quantitation of lipids was performed in the copresence of labeled and unlabeled lipids in a single DAG<sup>+</sup>-specific PIS spectrum. Briefly, the peak intensities of unlabeled lipid

molecular species were extracted and ratiometrically compared to the corresponding labeled ones as


$$C_{\text{labeled}}/C_{\text{unlabeled}} = I_{\text{labeled}}/I_{\text{unlabeled}},$$

$C_{\text{labeled}}$  and  $C_{\text{unlabeled}}$  represent the contents of individual species in the labeled and unlabeled lipid extracts, respectively, while  $I_{\text{labeled}}$  and  $I_{\text{unlabeled}}$  represent the sum of the heights of all the peaks of the isotopic distribution of the labeled ( $m/z$  of  $M - 1$ ,  $M$ ,  $M + 1$ ,  $M + 2$ ) and unlabeled ( $m/z$  of  $M$ ,  $M + 1$ ,  $M + 2$ ) lipid species, respectively. The results are presented as mean  $\pm$  SD. A <sup>13</sup>C correction, as described previously<sup>15</sup> with some modifications, was performed for quantification of phosphoinositides, to deal with the overlapping nature of the pseudomolecular ion peak of the species of interest with the  $M + 2$  isotope peak from another more saturated species that has a 2 Da lower mass (see Supporting Information Figure S2 for details).


**Safety Considerations.** The TMS-diazomethane was originally developed as a less toxic substitute for the compound diazomethane. However, inhalation of TMS-diazomethane may cause lung damage or central nervous system depression and should therefore be handled with the appropriate safety procedures in place, i.e., inside a fume hood, with adequate personal safety equipment. Excess TMS-diazomethane is to be neutralized using acetic acid. Care must be taken if handled in large volumes, as it is highly volatile and produces nitrogen gas.

## RESULTS AND DISCUSSION

**DAG<sup>+</sup>-Specific MPIS Enables Rapid and Simultaneous Screen of Phosphoinositide Molecular Species.** Clark et al. showed that TMS-diazomethane can be used to methylate the phosphate groups of phosphoinositides, thus providing a first solution to the issues of stability and sensitivity.<sup>11</sup> It has been shown previously that parts of the free hydroxyls in the inositol ring are methylated during the inositol phosphate methylation process. This produces phosphoinositides that are also methylated on the one or two free hydroxyls in the inositol ring (defined here as overmethylated phosphoinositide, in contrast to the completely methylated phosphoinositides that do not carry methyl groups on free hydroxyls in the inositol



**Figure 2.** Distribution of completed, overmethylated phosphoinositides in the DAG<sup>+</sup>-specific precursor ion scanning. (A and C) Simultaneous detection of completed, overmethylated 18:0/20:4-PtdInsP<sub>2</sub> and 18:0/20:4-PtdInsP<sub>3</sub> by precursor ion scanning at *m/z* of 627, which corresponds to the charged 18:0/20:4-DAG. (B and D) Statistical analysis showed significant similarity in the distribution of completed, overmethylated PtdInsP<sub>2</sub> or PtdInsP<sub>3</sub> of different fatty acyl chains, where the intensities of overmethylated PtdInsP<sub>2</sub> or PtdInsP<sub>3</sub> species with monomethylation of free hydroxyl groups in the inositol ring accounted for 50% ± 11.4% or 35% ± 1.5% of completed phosphoinositides, respectively.



**Figure 3.** Simultaneous detection of PtdInsP<sub>3</sub>, PtdInsP<sub>2</sub>, PtdInsP, and PtdIns, that have the same fatty acyl composition in a single DAG<sup>+</sup>-specific PIS spectrum. 18:0/20:4-PtdInsP<sub>3</sub>, PtdInsP<sub>2</sub>, PtdInsP, and PtdIns were simultaneously detected by precursor ion scanning at *m/z* of 627 that corresponds to the charged 18:0/20:4-DAG. A sequential mass interval of 108 Da can be found among PtdInsP<sub>3</sub>, PtdInsP<sub>2</sub>, PtdInsP, and PtdIns, which corresponds to the methylation of one phosphate. The derivatized PtdIns preferentially formed ammonium adducts.

ring). This may complicate the liquid chromatography based mass spectrometry analysis, due to the increasing complexity of phosphoinositide compositions. However, MS/MS analysis of completed or overmethylated phosphoinositide reveals that both yielded a fragment ion that corresponded to the charged diacylglycerol. For example, as shown in Figure 1, collision-induced dissociation (product ion analysis) of the completed or overmethylated 18:0/20:4-PtdInsP<sub>2</sub> (*m/z* of 1117 and 1131, respectively) can lead to the detection of the same fragment ion at *m/z* of 627, which corresponded to the charged 18:0/20:4-DAG. Presumably, the completed or overmethylated 17:0/20:4-PtdInsP<sub>3</sub> (*m/z* of 1211 and 1225, respectively) can yield the same daughter ion at *m/z* of 613, which corresponds to the charged 17:0/20:4-DAG. Therefore, phosphoinositides of various methylation states can be detected simultaneously by a precursor ion scan corresponding to the charged


diacylglycerol fragments, defined here as DAG<sup>+</sup>-specific PIS, by DIMS (Figure 2). More importantly, the profile of completed or overmethylated phosphoinositides in the precursor ion scan spectrum was of marked similarity among the phosphoinositides with different fatty acyl chains. For example, the intensity of overmethylated PtdInsP<sub>2</sub> and PtdInsP<sub>3</sub> species with monomethylation of free hydroxyl groups in the inositol ring were accounting for 50% ± 11.4% and 35% ± 1.5% of completed methylated ones, respectively (Figure 2). Similar to PtdInsP<sub>3</sub> and PtdInsP<sub>2</sub>, PtdInsP and PtdIns molecular species have their own unique profiling of completed or overmethylated phosphoinositides in the precursor ion scan spectrum (Supporting Information Figure S3).

Given that they have the same fatty acyl composition, the derivatized PtdInsP<sub>3</sub>, PtdInsP<sub>2</sub>, PtdInsP, or PtdIns yielded the



same DAG<sup>+</sup> fragment ion. This result indicates that PtdInsP<sub>3</sub>, PtdInsP<sub>2</sub>, PtdInsP, and PtdIns can be simultaneously and rapidly identified in a single DAG<sup>+</sup>-specific PIS. Figure 3 shows that 18:0/20:4-PtdInsP<sub>3</sub>, PtdInsP<sub>2</sub>, PtdInsP, and PtdIns was simultaneously presented with characteristic profiling of completed and overmethylated phosphoinositides in the DAG<sup>+</sup>-specific PIS spectrum. A sequential mass interval of 108 Da, corresponding to the methylation of one phosphate moiety, was also observed for 18:0/20:4-PtdInsP<sub>3</sub>, PtdInsP<sub>2</sub>, PtdInsP, and PtdIns. These results confirmed the feasibility of using DAG<sup>+</sup>-specific PIS for the rapid identification of particular phosphoinositide molecular species. Furthermore, the information mentioned above provides us higher confidence in the identification of low-abundant phosphoinositide species when compared to neutral loss scan analysis (described earlier by Clark et al.), since the peaks of low-abundant phosphoinositide species in NLS spectrum might be barely above the noise threshold. More importantly, it also could provide a significantly greater throughput compared to most existing methods.

In order to further test the feasibility of our novel method, a commercially available phosphoinositide mixture extracted from bovine brain was subjected to our DAG<sup>+</sup>-specific PIS analysis. DAG<sup>+</sup>-specific MPIS allowed rapid and sensitive identification of 32 PtdIns species, 28 PtdInsP species, 30 PtdInsP<sub>2</sub> species, and 3 PtdInsP<sub>3</sub> species from brain extracts. To our knowledge, this is the largest data set of phosphoinositides identified so far from mammalian brain extracts. These results also reveal the widest range of fatty acid composition in phosphoinositides ever found in mammalian brain extracts, although the precise reasons behind the need for such a wide variety of phosphoinositide species remain unknown. Previous studies provided a possible clue by showing that different phosphoinositide species could be produced upon stimulation using distinct ligands.<sup>9</sup> Therefore, our method, combined with specific agonist activation, allows us to begin the molecular dissection of various cellular outcomes regulated by differential molecular species of phosphoinositide. More importantly, the detection of three PIP<sub>3</sub> molecular species, 38:3PtdInsP<sub>3</sub>, 38:4 PtdInsP<sub>3</sub>, and 38:5 PtdInsP<sub>3</sub>, suggests that our DAG<sup>+</sup>-specific MPIS method is sufficiently sensitive to detect very low abundant phosphoinositide species (Supporting Information Figure S4A). In addition, many PtdInsP and PtdInsP<sub>2</sub> species not reported previously were also identified. For example, several PtdInsP/PtdInsP<sub>2</sub> species with long-chain fatty acyl, such as PtdInsP/PtdInsP<sub>2</sub> 42:5, PtdInsP/PtdInsP<sub>2</sub> 42:4, and PtdInsP<sub>2</sub> 42:2 (Supporting Information Figure S4B) were successfully detected. Furthermore, a large number of PtdInsP<sub>2</sub> and PtdInsP species containing highly polyunsaturated acyl chains were also found in our study, such as 36:8 PtdInsP<sub>2</sub>, 40:10 PtdInsP<sub>2</sub>, 40:9 PtdInsP<sub>2</sub> (Supporting Information Figure S4C). To confirm the identified species of phosphoinositides, the MS<sup>3</sup> analysis of specific DAG<sup>+</sup> ions was carried out to determine the fatty acyl composition of phosphoinositides. For example, as shown in Supporting Information Figure S5, both the fatty acyl chains of 40:10-PtdInsP<sub>2</sub> (*m/z* of 1133) were 20:5. Interestingly, we also found that the PtdIns, PtdInsP, and PtdInsP<sub>2</sub> species have similar fatty acid compositions (Figure 4). For example, the relative abundance (defined here as mol %) of PtdInsP/PtdInsP<sub>2</sub> species with identical fatty acyl composition within each class was similar. We speculate that this might provide an opportunity for cells to achieve rapid interconversions between




**Figure 4.** Relative abundance of different phosphoinositide molecular species identified from bovine brain extract within each class of PtdIns, PtdInsP, PtdInsP<sub>2</sub>, and PtdInsP<sub>3</sub>. The relative abundance of each phosphoinositide species, defined here as mol %, was calculated by dividing its peak intensity, including the overmethylated ones, by the summing peak intensity of phosphoinositide species detected within each class.


different phosphoinositide classes, though whether such rapid interconversions are really existed is still unclear. However, to determine whether such similarity is caused by interconversions occurring as a function of sample handling or in vivo processing, we designed experiment as shown in Supporting Information Figure S6. As shown in Supporting Information Figure S7, no conversions of phosphoinositide standards to more or less phosphorylated structures were detected in the three group samples. These results suggest that the interconversions occurring as a function of sample handling is neglectable in our case.

**Profiling of Phosphoinositides from Different Cancer Cell Extracts.** We next tested whether the DAG<sup>+</sup>-specific MPIS mentioned above can be used to analyze phosphoinositides directly extracted from cell cultures. To this end, three prostatic cancer cell lines, namely, PC3, C4-2B, as well as Du145, were subjected to a rapid, yet comprehensive screening of phosphoinositides using our DAG<sup>+</sup>-specific MPIS. As a result, we identified ~90 phosphatidylinositol and phosphatidylinositol phosphate molecular species (Supporting Information Table S2). Supporting Information Figure S8 showed that the range of fatty acid chains in phosphoinositides of these cancer extracts was similar to that in brain extracts. However, the relative abundance of different phosphoinositide species in cultured cells differed from those analyzed from brain extracts. For example, the most abundant species detected in bovine brain extract were 38:4 PtdInsP<sub>2</sub>, 38:3 PtdInsP<sub>2</sub>, 34:1 PtdInsP<sub>2</sub>, and 36:1 PtdInsP<sub>2</sub>, comprising ~37, ~17, ~8, and ~6 mol % of all identified PtdInsP<sub>2</sub> species, respectively (Supporting Information Table S3). In contrast, the most abundant phosphoinositide species detected in C4-2B cell line were 38:4, 36:1, 36:2 and 38:3 fatty acyl species, which account for ~19, ~17, ~14, and ~12 mol % of all PtdInsP<sub>2</sub> species identified. Bovine brain has conspicuously higher levels of 38:4 series phosphoinositides compared with C4-2B cells, possibly because the serum used in cell culture is relatively low on 20:4 (arachidonate) or its metabolic precursor 18:2 (linoleic acid). This result is consistent with previous observation.<sup>11</sup>

Taken together, DAG<sup>+</sup>-specific MPIS provide a simple, fast, sensitive, and high-throughput approach for comprehensive



**Figure 5.** Scheme for isotopic labeling of phosphoinositides by phosphate methylation based on acid-catalyzed H/D exchange and methanolysis of TMS-diazomethane.




**Figure 6.** Efficiency of isotopic labeling of phosphoinositides with phosphate methylation by diazomethane generated in situ. No unlabeled phosphoinositides were detected in the isotopically labeled samples, demonstrating the labeling was close to complete. \* indicates skewed isotope distribution with one minor peak of a mass difference of 1 Da.

profiling of phosphoinositide species from different biological sample, which should accelerate the elucidation of further molecular species of phosphoinositides across different cells, tissues and organisms.

**Isotopic Labeling of Phosphoinositides Based on Phosphate Methylation.** Previously, Kuhnel et al. showed that methyl-esterification of carboxylic acids by TMS-diazomethane/MeOH proceeds through in situ methanolytic liberation of diazomethane.<sup>16</sup> In this reaction, carboxylic acid acts as a catalyst for the methanolysis of TMS-diazomethane, first as a general acid, then as a general base, in the generation of free diazomethane. It was shown that, by isotopic labeling, deuterated diazomethane (CD<sub>2</sub>N<sub>2</sub>) can be generated by acid-catalyzed H/D exchange and methanolysis, with the MeOD acting as a deuterium reservoir. The CD<sub>2</sub>N<sub>2</sub> can then be used to label the carboxylic acid with -CHD<sub>2</sub> group. In addition, TMS-diazomethane enables relatively fast and clean esterification of protonated phosphate groups of phosphoinositides. We therefore reasoned that the phosphate methylation by TMS-diazomethane may behave in a similar way to that of carboxylic acid (Figure 5). These findings encouraged us to test MeOD/DCl as a D reservoir during the isotopic labeling of phosphoinositides by TMS-diazomethane. It is reasonable to assume that when D-labeled heavy diazomethane was utilized instead of a light reagent, a mass increase should be detected by the mass analyzer, with such an increase being directly proportional to the amount of available phosphate groups in

the phosphoinositides. Therefore, the larger the number of phosphate groups in a phosphoinositide, the larger the mass difference between the <sup>1</sup>H- and the D-labeled isotopic pair can be identified. As shown in Figure 6, rapid and efficient phosphate methylation and isotopic labeling of phosphoinositides, either PtdInsP, or PtdInsP<sub>2</sub>, or PtdInsP<sub>3</sub>, with -CHD<sub>2</sub> groups was achieved. The mass differences between the <sup>1</sup>H- and the D-labeled isotopic pairs of PtdInsP, PtdInsP<sub>2</sub>, and PtdInsP<sub>3</sub> were 6, 10, and 14 Da, respectively. No unlabeled phosphoinositides were detected in the isotopically labeled samples, demonstrating that the process of generating CD<sub>2</sub>N<sub>2</sub> and phosphate methylation of phosphoinositides was achieved close to completion. However, it should be noted that the isotopic cluster of D-methylated phosphoinositide shows the presence of one less abundant peak that is 1 mass unit smaller than the monoisotopic peak for the D-methylated phosphoinositide. For example, the labeled 38:4 PtdInsP and PtdInsP<sub>2</sub> (m/z of 1015 and 1127, respectively) exhibited a skewed isotope distribution and one minor peak with mass difference of 1 Da (m/z of 1014 and 1126, respectively) (Figure 6, indicated by an asterisk). This is partly due to a small percentage of <sup>1</sup>H atoms in the D-marked MeOD reagent used in our experiments, as well as due to the effects of back-exchange during the generation of D-diazomethane. In addition, a small proportion of the resulting methyl ester would be -CD<sub>3</sub> group, due to the H/D exchange of acidic protons from the phosphate group with deuterium atoms that in the environment, because



**Figure 7.** Precursor ion scan spectra showing the isotopic pairs of H- and D-labeled phosphoinositides at known ratios. The precursor ion scan for the 18:0/20:4-PtdInsP<sub>2</sub> from phosphoinositides mixing ratios 1:1, 2:3, 1:3, 3:1, 1:10, and 10:1 are shown in panels A–F. The monoisotopic *m/z* for the H- and D-labeled 18:0/20:4-PtdInsP<sub>2</sub> are 1117 and 1127, respectively. The H/D ratios were calculated from the sum of the heights of all the peaks of the isotopic distribution of each component of the isotopic pair. The figure illustrates the ease of identification of matched H- and D-labeled pairs and the option for quantitation using the ratio of phosphoinositide intensities from the precursor ion scan spectra.

previous studies had shown that two of the protons in the resulting methyl ester originate from the diazomethane and the other one is the “donated” acidic proton from the carboxylic acid (phosphate, in our case).<sup>16</sup> Nonetheless, our results show that the use of TMS-diazomethane, with MeOD/DCl as a D reservoir, is feasible for isotopic labeling of phosphoinositides, which could then be employed for relative quantitation of phosphoinositides.


#### Development and Validation of the Quantitation Method.

Quantitation of phosphoinositides using our method was simplified by the copresence of labeled and unlabeled phosphoinositides in a single DAG<sup>+</sup>-specific PIS spectrum, as the “light” and “heavy” labeled phosphoinositides yielded an identical DAG<sup>+</sup> fragment ion during collision-induced dissociation (Supporting Information Figure S9). Therefore, it is possible to obtain ratios of phosphoinositides isolated at different conditions directly from the DAG<sup>+</sup>-specific PIS spectrum, as differences in the ion abundances between the isotopic pairs should reflect differences in the concentration between the two conditions. The ratios can then be calculated from the peak heights of the monoisotopic masses of each component of the isotopic pair, expressed here as H/D ratios. However, more cautions should be taken in our case. First, whether the overmethylated phosphoinositides to be included for the quantitative analysis? Our results show that, in most cases, the distribution of deuterium-labeled completed or overmethylated PtdInsP, PtdInsP<sub>2</sub>, and PtdInsP<sub>3</sub> is similar to that of hydrogen-labeled ones. For instance, the Supporting Information Figure S10 shows that ratios of hydrogen/deuterium-labeled overmethylated phosphoinositides with different fatty acyl chains were parallel with that of completed-methylated ones. In addition, as shown in Supporting Information Figure S10, there was an overlap between the <sup>1</sup>H-labeled PtdInsP<sub>2</sub> with two overmethylations and D-labeled PtdInsP<sub>2</sub> with one overmethylation, which might turn the quantitative analysis into a more complicated procedure. Furthermore, the overmethylated phosphoinositide

peaks are more prone to noise interference than that of completed-methylated ones, especially for low-abundant phosphoinositides. Meanwhile, due to the paired comparison nature of isotopic labeling based relative quantitative analysis, overmethylated phosphoinositide peaks can be easily included for quantitative analysis if necessary. Therefore, we preferred to use the completed-methylated phosphoinositides for quantitative analysis. Second, as mentioned above, the height of the monoisotopic peak in the D-labeled phosphoinositide was affected by the presence of a small proportion of <sup>1</sup>H in the isotopic labeling reagent, as well as the effect of back-exchange during the generation of CD<sub>2</sub>N<sub>2</sub> and H/D exchange of acidic proton from the phosphate group. An alternative way to calculate the H/D ratios was then taken, by calculating the sum of the heights of all the peaks of the isotopic distribution of each component of the isotopic pair. For example, the sum intensity of the peaks corresponding to unlabeled (*m/z* of 1117, 1118 and 1119) and labeled (*m/z* of 1126, 1127, 1128 and 1129) 18:0/20:4-PtdInsP<sub>2</sub> was taken to calculate the H/D ratio for unlabeled and labeled 18:0/20:4-PtdInsP<sub>2</sub> (monoisotopic *m/z* of 1117 and 1127, respectively). In addition, the derivatized PtdInsP<sub>3</sub>, PtdInsP<sub>2</sub>, PtdInsP, or PtdIns that has the same fatty acyl composition, regardless of the labeled or unlabeled forms, yielded a same DAG<sup>+</sup> fragment ion. This allowed the simultaneous quantification of PtdInsP<sub>3</sub>, PtdInsP<sub>2</sub>, PtdInsP, and PtdIns in a single precursor ion scan (Supporting Information Figure S11).

To characterize the analytical properties and the accuracy of quantitation by phosphate methylation with CD<sub>2</sub>N<sub>2</sub>, the PC3 cell cultures were divided into several aliquots (~5 × 10<sup>5</sup> cells/aliquot) and separated into following proportions of cell volume, 10:1, 3:1, 3:2, 1:1, 1:3, and 1:10, before phosphoinositide extraction and separate labeling with either <sup>1</sup>H- or D-labeled diazomethane. Our results showed that the observed H/D ratios were consistent with the theoretical ratios of PtdInsP<sub>3</sub>, PtdInsP<sub>2</sub>, and PtdInsP present in the mixture (Supporting Information Table S4); the average CV from

three replicates was 18.9%. Figure 7 shows that the difference in signal intensity of the isotopic pair of the 38:4 PtdInsP<sub>2</sub> was found to be in good agreement with the theoretically expected ratios, with an average CV of less than 12% from three replicates. To further illustrate the linearity of the actual ratios, measurements corresponding to data for the phosphoinositides at different H/D proportions were plotted against their respective theoretical values. The resulting graph (Figure 8)



**Figure 8.** Comparison of the measured H/D ratios with their respective theoretical values. A straight line ( $r^2 = 0.99$ ) was observed when the H/D ratios of the PtdInsP<sub>2</sub> from PC3 cell lines at different H/D ratios were graphed against their theoretical H/D mix ratios. Error bars represent the difference between the highest and the lowest values measured for each ratio.

exhibited linearity with an  $R^2$  of 0.99 and a slope of 1.03 for PtdInsP<sub>2</sub> species over a dynamic range that was close to 2 orders of magnitude. Similar linearity was observed for the isotopic pairs of the PtdInsP (Supporting Information Figure S12). For quantification of PtdInsP<sub>3</sub>, caution is required. As shown in Figure 5, there would be an overlap between the D-labeled PtdInsP<sub>3</sub> (+14D, due to 14 H/D exchange happened in seven methyl groups) and the unlabeled but overmethylated PtdInsP<sub>3</sub> (+CH<sub>2</sub>) in the DAG<sup>+</sup>-specific PIS spectrum. For example, the exact  $m/z$  value of the labeled 38:4 PtdInsP<sub>3</sub> and unlabeled overmethylated 38:4 PtdInsP<sub>3</sub> is 1239.65 and 1239.58, respectively. Therefore, the raw ratios of PtdInsP<sub>3</sub> were corrected with a factor of  $-0.35$ , as the intensity of the overmethylated form accounting for  $\sim 35\%$  of completed methylated PtdInsP<sub>3</sub> (Figure 2). Alternatively, the quantification of PtdInsP<sub>3</sub> can be carried out by high-resolution mass spectrometry, which is of the power to resolve the tiny mass difference mentioned above.

**Potential Pitfalls in Quantitation of Phosphoinositides.** Although very high accuracy for quantitation of a given phosphoinositides is possible under optimal conditions, there are several situations where quantitation accuracy could be compromised. An important issue that is common to all MS-based methods for quantitation is the signal-to-noise ratio (S/N) for each of the two peaks used for comparison. For example, low-abundance phosphoinositide peaks may be barely above the noise threshold, which may receive an inordinate contribution of signal intensity from background “noise”; the apparent ratio between the two peaks will therefore be less or more than the actual ratio. In contrast to metabolic labeling methods, efficiency of labeling, sample preparation biases introduced by slightly different purification conditions, such as

differential losses of certain phosphoinositide subclasses in the organic phase during their extraction, may affect the phosphoinositide ratios, because the samples are mixed after separately extraction and labeling. However, such quantitation error can be corrected by inclusion of one ISD prior to lipid extraction and labeling. In addition, a strategy by considering both forward and reverse labeling for the control and treated samples for all the duplicates, as shown in Supporting Information Figure S13, was suggested for the relative quantitative analysis. This analysis could neglect effects coming from the usage of different labels.

## CONCLUSIONS

The results presented here demonstrate the feasibility of rapid profiling and relative quantitation of phosphoinositides by DAG<sup>+</sup>-specific MPIS based on phosphate methylation and isotopic labeling by chemical derivatization with diazomethane. Our novel assay is not only rapid, simple, and sensitive, but also offers important advantages over existing methods for phosphoinositide measurement. First, the DAG<sup>+</sup>-specific MPIS enables the simultaneous identification of PtdInsP<sub>3</sub>, PtdInsP<sub>2</sub>, PtdInsP, and PtdIns in an efficient and time-saving manner, which provide a significantly greater throughput compared to most existing methods. Using our DAG<sup>+</sup>-specific MPIS strategy, we were able to identify more than 90 phosphoinositides from different biological samples, such as bovine brain extract as well as prostatic cancer cell cultures. Our DAG<sup>+</sup>-specific MPIS method, combined with isotopic labeling by diazomethane allows for simultaneous and precise relative quantification of PtdInsP<sub>3</sub>, PtdInsP<sub>2</sub>, PtdInsP, and PtdIns, without the need of multiple ISDs. Isotopic labeling based quantitation of phosphoinositides also minimizes matrix effects or run-to-run differences that could potentially affect results. The remaining limitation to identify regioisomers of polyphosphoinositides could be resolved by combining our method with the appropriate chromatographic separation methods.

## ASSOCIATED CONTENT

### Supporting Information

Additional information as noted in text. This material is available free of charge via the Internet at <http://pubs.acs.org>.

## AUTHOR INFORMATION

### Corresponding Author

\*E-mail: [fqyang@ibp.ac.cn](mailto:fqyang@ibp.ac.cn). Phone: +86-10-6488-8581. Fax: +86-10-6488-8581.

### Author Contributions

<sup>†</sup>T.C. and Q.S. contributed equally to this work.

### Notes

The authors declare no competing financial interest.

## ACKNOWLEDGMENTS

We thank all members of the Yang lab for critical reading of the manuscript and excellent technical support. We thank Dr. Torsten Juelich for excellent language editing. This work was supported by grants from the National Basic Research Program of China (Grant Nos. 2010CB833703, 2012CB966803, 2011CB915501, and 2014CBA02003), the National Natural Science Foundation of China (Grant Nos. 90919047, 81028009, and 31100614), and Novo Nordisk-CAS Research Foundation NNCAS-2011-1.



## ■ REFERENCES

- (1) Martin, T. F. *Annu. Rev. Cell Dev. Biol.* **1998**, *14*, 231–264.
- (2) Di Paolo, G.; De Camilli, P. *Nature* **2006**, *443*, 651–657.
- (3) Jones, D. R.; Varela-Nieto, I. *Mol. Med.* **1999**, *5*, 505–514.
- (4) Pendaries, C.; Tronchere, H.; Plantavid, M.; Payrastre, B. *FEBS Lett.* **2003**, *546*, 25–31.
- (5) Balla, T. *Physiol. Rev.* **2013**, *93*, 1019–1137.
- (6) Guillou, H.; Stephens, L. R.; Hawkins, P. T. *Methods Enzymol.* **2007**, *434*, 117–130.
- (7) Jones, D. R.; Ramirez, I. B.; Lowe, M.; Divecha, N. *Nat. Protoc.* **2013**, *8*, 1058–1072.
- (8) Wenk, M. R.; Lucast, L.; Di Paolo, G.; Romanelli, A. J.; Suchy, S. F.; Nussbaum, R. L.; Cline, G. W.; Shulman, G. I.; McMurray, W.; De Camilli, P. *Nat. Biotechnol.* **2003**, *21*, 813–817.
- (9) Milne, S. B.; Ivanova, P. T.; DeCamp, D.; Hsueh, R. C.; Brown, H. A. *J. Lipid Res.* **2005**, *46*, 1796–1802.
- (10) Ivanova, P. T.; Milne, S. B.; Myers, D. S.; Brown, H. A. *Curr. Opin. Chem. Biol.* **2009**, *13*, 526–531.
- (11) Clark, J.; Anderson, K. E.; Juvin, V.; Smith, T. S.; Karpe, F.; Wakelam, M. J.; Stephens, L. R.; Hawkins, P. T. *Nat. Methods* **2011**, *8*, 267–272.
- (12) Ong, S. E.; Blagoev, B.; Kratchmarova, I.; Kristensen, D. B.; Steen, H.; Pandey, A.; Mann, M. *Mol. Cell. Proteomics* **2002**, *1*, 376–386.
- (13) Ross, P. L.; Huang, Y. N.; Marchese, J. N.; Williamson, B.; Parker, K.; Hattan, S.; Khainovski, N.; Pillai, S.; Dey, S.; Daniels, S.; Purkayastha, S.; Juhasz, P.; Martin, S.; Bartlett-Jones, M.; He, F.; Jacobson, A.; Pappin, D. J. *Mol. Cell. Proteomics* **2004**, *3*, 1154–1169.
- (14) Gray, A.; Olsson, H.; Batty, I. H.; Priganica, L.; Peter Downes, C. *Anal. Biochem.* **2003**, *313*, 234–245.
- (15) Han, X.; Gross, R. W. *Mass Spectrom. Rev.* **2005**, *24*, 367–412.
- (16) Kuhnle, E.; Laffan, D. D.; Lloyd-Jones, G. C.; Martinez Del Campo, T.; Shepperson, I. R.; Slaughter, J. L. *Angew. Chem., Int. Ed.* **2007**, *46*, 7075–7078.